

Price List G 2008
Diagnostic Microscope
MÖLLER SPECTRA
with Carrying Systems
FS/CU/WU 1-11

Part no.	Description	Weight approx. kg	Qty.	Price €
655 438	<p>Stereo Operating Microscope MÖLLER SPECTRA 300 - Dental</p> <p>-----</p> <p>Microscope in swivable and tilt-able single-arm suspension; manual 5-step magnification, T-handle, 3 step slide incl. yellow filter, 45 degrees eyepiece head, convergent with 2 wide angle oculars, 12,5 x with fixation; front lens f = 250 mm. Operating Manual</p>		1	
655 436	<p>Stereo Examination Microscope MÖLLER SPECTRA 500 - ENT</p> <p>-----</p> <p>Microscope in swivable and tilt-able single-arm suspension; manual 5-step magnification; arch handle, 5 step slide incl. green filter, 0 degrees eyepiece head, convergent with 2 wide angle oculars 12,5x with fixation; front lens f= 250 mm. Operating Manual</p>		1	
655 439	<p>Stereo Operating Microscope MÖLLER SPECTRA 900 - Ophthalmic</p> <p>-----</p> <p>Microscope in swivable and tilt-able single-arm suspension; manual 5-step magnification; straight handles; 5 step slide incl. yellow filter, 45 degr. convergent eyepiece head with 2 wide angle oculars 12,5x with fixation; front lens f = 200 mm. Operating Manual</p>		1	
<p><u>CARRYING SYSTEMS</u></p>				
615 480	<p>Ceiling Unit CU 1-11/150 W</p> <p>-----</p> <p>Counter weight balanced spring joint arm with rotatable carrier. Length: 1.480 mm, Lift 600 mm; Light fibre cable; Continuously adjustable cold light illumination (150 W Halogen illumination with 2 bulbs and quick switch over) Power supply: 110/230 V.</p>		1	
604 555	<p>Ring for fixation of Ceiling Unit</p>		1	

Part no.	Description	Weight approx. kg	Qty.	Price €
615 460	<p>WALL UNIT -----</p> <p>Wall mount WU 1-11 / 150 W -----</p> <p>Counter weight balanced spring joint arm with rotatable carrier, length: 1.370 mm, lift 380 mm; light fibre cable; Continuously adjustable cold light illumination (150 W Halogen illumination with 2 bulbs and quick switch over) Power supply: 110/230 V incl. - power cable, length 4 meters, with German plug, - Wall plate P/N 615 461</p>		1	
615 433	<p>CONNECTING CABLES -----</p> <p>Power cable for WU 1-11, length 4 m, German plug</p>		1	
615 440	<p>FLOORSTAND -----</p> <p>Floorstand FS 1-11 for MÖLLER SPECTRA and DENTA 300 -----</p> <p>Stand with 4 rollers; weight counter balanced spring joint arm with rotatable carrier. Continuously adjustable cold light illumination (150 W Halogen illumination with 2 bulbs and quick switch over).</p>		1	
614 215	<p>CONNECTING CABLES for FS 1-1x -----</p> <p>Power cable 6 m with German plug</p>		1	
614 220	<p>Power cable 6 m with US plug</p>		1	

Part no.	Description	Weight approx. kg	Qty.	Price €
	ACCESSORIES: For MÖLLER SPECTRA			
655 325	Front Lens f = 175 mm		1	
655 320	Front Lens f = 200 mm		1	
655 327	Front Lens f = 250 mm		1	
655 328	Front Lens f = 300 mm		1	
655 329	Front lens f = 350 mm		1	
655 330	Front Lens f = 400 mm		1	
655 334	Front lens revolver for quick change of 2 front lenses (without front lenses)		1	
655 393	Focusing device f=200 mm for SPECTRA, incl. front lens		1	
655 394	Fine focusing device f = 250 mm for SPECTRA, incl. front lens		1	
655 398	Fine focusing device f = 400 mm for SPECTRA, incl. front lens		1	
655 399	Focusing device f=300 mm for SPECTRA, incl. front lens		1	
655 303	Filter disc with green filter		1	
655 348	Filter slide with orange filter for SPECTRA.		1	
655 300	Blue filter for SPECTRA		1	
655 310	Eyepiece Head, 0°, convergent		1	
655 335	0° Eyepiece Head, convergent, with PD spindle, 12,5 x for SPECTRA		1	
655 312	Eyepiece Head, 45°, convergent		1	
655 337	45° Eyepiece Head, convergent, with PD-spindle, 12,5 x		1	
655 287	Handle (double arche)		1	
655 301	Handle (double straight) for SPECTRA 900		1	
655 391	Adapter for sterile microscope cover P/N 190 304 - SPECTRA/DENTA		1	

Part no.	Description	Weight approx. kg	Qty.	Price €
190 304	Sterile Microscope Cover, *) type ST, length 140 cm (min. quantity to order = 10 pc.) Charge-No.:		1	
655 341	Attachment for 1-chip lipstick TV-camera (mount M 15,5 x 0,5)/ SPECTRA		1	
655 346	TV attachment for TV camera 1/2" (c-mount)/SPECTRA		1	
	SPARE PARTS: for FS 1-11 =====			
614 215	Power cable 6 m with German plug		1	
628 115	Cap 12 D x 80, sterilizable		1	
123 112 61	Spare lamp for cold light source Halog-Lp-Refl-G-15V-150W for carrying systems 1100, 2012, 3-21/3-31 and 1-1x		1	
126 401 07	Light fibre cable for *) Carrying Systems FS 1100, 2011, 2012, 3-21, length 2,15 m		1	
126 401 54	Light fibre cable for *) FS/WU/CU 1-1x, CU 1100 / WU 1100 Length 2,45 m		1	
	Constructional work on site: ----- - Ceiling Unit: For further installation details please refer to our "Technical Information Ceiling Unit" - Wall Mount: see our "Technical Information" A plug socket with the requested voltage must be available. - Floorstand: A plug socket with the requested voltage must be available. *) = OEM parts Prices are to be understood net ex works Wedel, incl. packing. Validity until end of 2008. Constructional modifications and price changes are reserved.			